

# Densiron® XTRA

Heavy silicone oil  
with molecularly design

New  
Compo-  
sition



Unique molecularly design

Heavier than water

Highly resistant against emulsification<sup>1</sup>

Ideal for inferior pathologies

25G compatible

1 Caramoy A, Schröder S, Fauser S, Kirchhof B (2010) In vitro emulsification assessment of new silicone oils. Br J Ophthalmol 94, 509-512

**Syringe** G-80925 Densiron® Xtra 10 ml  
Syringe, 1 pc. per box, sterile

Fluoron GmbH reserves the right to make changes to compositions in response to recent developments. Fluoron GmbH and GEUDER AG do not assume liability for the accuracy of given statements. Regarding product availability in the single countries please contact your distributor or GEUDER AG.

**Manufacturer:** Fluoron GmbH Magirus-Deutz-Strasse 10 89077 Ulm Germany  
Phone: +49 731 205 5997 0 Fax: + 49 731 205 5997 28 info@fluoron.de www.fluoron.de

**Distributed by:** GEUDER AG Hertzstrasse 4 69126 Heidelberg Germany  
Phone: +49 6221 3066 Fax: +49 6221 303122 info@geuder.de www.geuder.de

**Geuder®**  
Precision made in Germany

**FLUORON®**  
Leading in purity and variety