

## VERRE GRAND CHAMP À USAGE UNIQUE SINGLE USE WIDE FIELD LENS

FDA APPROVED

### "Panoramic"



#### **Application Principale:**

Idéal pour les décollements de rétine et les déchirures rétiniennees géantes.

#### **Caractéristiques et Avantages :**

Petite hauteur facilitant la manipulation sur l'orbite.  
Disponible avec un anneau silicone et manche, amovible.  
Grand champ de vision de l'ensemble de la rétine.  
Utilisable avec un inverseur à image.  
Boite de 5 stériles

#### **Main application:**

Ideal for retinal detachments and giant retinal tears.

#### **Feature and advantage:**

Small height for easy handling on the orbit.  
Available with removable handle and silicone ring.  
Wide field of vision of the entire retina.  
Used with an image inverter.  
Box of 5 sterile



MV LV48 Lentille - Lens	Champ de Vision Field of View	Grossissement Magnification
Verre Grand Champ à usage unique Wide field lens single use	130 / 150°	0,39x

## PHAKOS

62 rue Kléber 93100 MONTREUIL-FRANCE  
Tél. +33 (0)1 49 88 74 00 - Fax.+33 (0)1 49 88 72 88  
phakos@phakos.com

CE 0459 - Fabriqué en France