

vitraTM probes


Premium Next Generation Laser Probes



20G - 23G - 25G

vitraTM probes

Premium Next Generation Laser Probes



Complete range

- 20, 23, 25 Gauge straight probes
- 20 Gauge curved probe
- 23, 25 Gauge flexible curved probes
- 23, 25 Gauge steerable probes



Comfort of use

- Optimized handpieces for an ideal probes handover
- Super flexible fiber cords
- Color coded by gauge size



23G and 25G Steerable Probes

- Steerable concept: pressure-controlled curvature from the handpiece
- 0-90° adjustable curve
- Unparalleled access to peripheral retina
- 360° rotatable handle
- Trocar cannula compatible



23G and 25G Flexible Curved Probes

- Flexible tip with memory shape effect
- Exclusive design for an optimal use through 23G and 25G trocars

Compatible with a large range of laser's connectors

TYPE 1	TYPE 2	TYPE 3
Quantel Medical™ <i>Vitra/Supra</i>	Quantel Medical™ <i>Viridis</i>	Iridex® Bausch & Lomb®
	Alcon®/Coherent®/Zeiss®/Ellex®	

©2015. Quantel Medical, Vitra Probes is a trademark of Quantel Medical. All rights reserved.

www.quantel-medical.com

Headquarters

Quantel Medical
11, rue du Bois Joli - CS40015
63808 Cournon d'Auvergne - FRANCE
Tel: +33 (0)4 73 745 745
Fax: +33 (0)4 73 745 700
E-mail: contact@quantel-medical.fr

North America

Quantel USA
601 Haggerty Lane
Bozeman, MT 59715 – USA
Tel: +1 877 782 6835
Fax: +1 406 522 2005
E-mail: info@quantelmedical.com

Representative Offices

Thailand, Chiang Mai
Brazil, Rio De Janeiro

 Quantel
medical

CE
0050