

PRISME TONOMÈTRE À USAGE UNIQUE SINGLE USE TONOMETER PRISM

"Tonoclean"



Sécurité :

Élimine le risque d'infection croisée entre les patients et les protège des maladies transmises par l'intermédiaire du film lacrymal.

Facilité d'utilisation :

Les prismes Tonoclean sont stériles emballés individuellement, permettant le montage sans contact direct avec le prisme.

Type :

Jetables, stériles, emballés individuellement
Montage : convient au tonomètre type Goldman
Longueur : 30mm
Emballé : Boite de 100

Safe

Eliminates the risk of cross-infection between patients and protects them against crossinfection of diseases transmitted through the tear film.

Ease of use

Tonoclean prisms are sterile, individually wrapped, allowing mounting without direct contact with the prism.


Type:

Disposable, sterile, individually wrapped
Mounting : suitable for Goldman tonometer type
Length : 30mm
Packaged : Box of 100

MV CT50 Prisme Tonomètre – Tonometer Prism

PHAKOS

62 rue Kléber 93100 MONTREUIL-FRANCE
Tél. +33 (0)1 49 88 74 00 - Fax.+33 (0)1 49 88 72 88
phakos@phakos.com

 CE 0459 - Fabriqué en France